

Micro

MicroRint

Medical

a subsidiary of  
VIASYS Healthcare

A simple method to  
measure airway resistance



## Medical

## A simple test for childhood asthma

Micro Medical's unique design has taken a hitherto specialised measurement out of the pulmonary function laboratory and into the clinic or home. MicroRint enables airway resistance to be measured with the same ease as peak flow but without effort or technique from the patient. The subject simply breathes passively through a mouthpiece or face mask. A rapidly occluding valve automatically interrupts the airflow for a period so brief as to be imperceptible to the patient. The airway resistance result is automatically computed and displayed.

The whole procedure takes only a few minutes and can be performed on all age ranges from neonates to adults. With the transducer weighing only 350g, the instrument is palm sized, easy to use and comes complete with all accessories in a sturdy carrying case.



### MicroRint Cat. No. MR5000

#### Features

- Lightweight and portable, operating from internal rechargeable batteries.
- Rapid, non-invasive inspiratory and expiratory measurements.
- An alternative to spirometry for pre-school age children.
- Well tolerated by neonates and geriatric patients.
- Fully configurable features menu.
- Pre and post medication comparisons with predicted values.
- Results printed directly to Hewlett Packard or Canon printers. Compatible models specified by Micro Medical.
- Rida airway resistance database and analysis software available as an option.
- 1000 + patient memory.
- The MicroRint transducer (Cat. No. MRT6000) can be purchased as an option for SuperSpiro, MicroLoop and MicroLab spirometers.




## Rint database and analysis software

An option to MicroRint is Rida, a 32 bit PC based software for displaying results and the waveforms developed during airway resistance testing

Rida has a user friendly, modern, multi-window visual interface and is compatible with all the latest Windows operating systems. Patient databases are easily created, a powerful search facility and the long term trending of patient results are also possible.

Bronchodilator and Steroid assessment can be made and Flow/Volume, Pressure/Time or Flow/Time waveforms can be displayed simultaneously ensuring quality testing and accurate results are obtained different sets of normal values are selectable and the printout format is fully customisable



### System Requirements

- PC with Microsoft Windows 2000 or later
- 32 MB of RAM
- 50 MB of free hard drive space
- Free Serial Port or USB Port
- Microsoft® Internet Explorer version 4.01 or above
- Video: SVGA 800x600, 256 colour

Rida Cat. No. RD1000


### Features

- Multi-window layout for ease of use
- Real-time waveform displays ensuring quality testing
- Bronchodilator and Steroid assessment is possible
- Choice of predicted values
- Configurable print out format



## The problem

Wheeziness is a common problem in pre-school children. Objective assessment of this is difficult because such children cannot usually co-operate with conventional tests of pulmonary function such as peak expiratory flow (PEF), or spirometry (FEV<sub>1</sub> - FVC). This has hindered both the diagnosis of respiratory tract conditions and assessment of the response to therapy.



**Q** Respiratory flow rate just before interruption

**PaO** Airway pressure (mouth pressure) immediately after interruption

**Pa** Alveolar pressure. Equilibrates with PaO immediately after interruption

**Rint** Airway resistance result (kPa/l/s) computed from PaO/Q

## Specification

|                       |  |
|-----------------------|--|
| Storage Capacity      | 1024 tests   |
| Printer compatibility | Hewlett Packard or Canon printer<br>[adaptor needed - Cat. No. MLA 350]                                    |
| Predicted values      | Percent predicted and normal range reported  |
| Display               | Graphic LCD 240 x 160 pixels   |
| Transducer type       | Flow: Pneumotachometer   |
| Pressure              | Piezo resistive  |
| Resolution            | 0.01 l/s   |
| Flow range            | 0 to 2 l/s   |
| Power supply          | Either primary 230V AC 50 Hertz secondary 9V DC 300mA<br>or primary 120V AC 60 Hertz secondary 9V DC 300mA |
| Battery Pack          | Micro controller unit: Rechargeable NiCad 3.6V 600mA hours<br>Transducer: Rechargeable NiCad 6V 50mA hours |
| Dimensions            | Micro controller unit 274 x 134 x 36mm<br>Transducer: 160 x 65 x 40mm                                      |
| Weight                | 1kg, 2.75kg with carrycase and accessories   |

## Bibliography

- Isobel Dundas and Sheila McKenzie: Measurement of lung function in preschool children *Airways J* 2003; 1: 30-3
- Bridge PD, Ranganathan S, McKenzie SA: Measurement of airway resistance using the interrupter technique in pre-school children in the ambulatory setting. *European Respiratory Journal*, 1999; 13, 1-5
- Kannisto S, Vanninen E, Remes K and Korppi M: Interrupter technique for evaluation of exercise-induced bronchospasm in children. *Pediatric Pulmonology*, vol. 27, Number 3, March 1999; 203-207
- Phagoo SB, Wilson NM, Silverman M: Evaluation of a new interrupter device for measuring bronchial responsiveness and the response to bronchodilator in 3 year old children. *European Respiratory Journal*, 1996; 9, 1374-1380
- Bridge PD, Lee H, M Silverman: A portable device based on the interrupter technique to measure bronchodilator response in schoolchildren. *European Respiratory Journal*, 1996; 9, 1368-1373
- Chowienzyk PJ, Lawson CP, Lane S et al. A flow interruption device for the measurement of airway resistance. *European Respiratory Journal*, 1991; 4, 623-628
- Enrico Lombardi MD, Elio Novembre, MD and Alberto Vierucci MD: Pediatric Allergy and Pulmonology Center, Anna Meyer Children's Hospital, University of Florence, Italy

The MicroRint and Rida are part of an extensive range of respiratory monitoring equipment manufactured by MicroMedical Ltd and is offered as Cat No. MR5000. (Rida Cat No. RD1000).

Micro Medical Ltd pursues a policy of continuing improvement in design, production and performance of its products; the right is therefore reserved to vary details at any time and without notice.



### Micro Medical Limited

PO BOX 6, Rochester, Kent, ME1 2AZ, UK

Telephone 01634 893500

Fax 01634 893600

International +44 1634 893500

Email [micromedical@viasyshc.com](mailto:micromedical@viasyshc.com)

[www.micromedical.co.uk](http://www.micromedical.co.uk)



Local Representative:

