Sys*Stim 294 Specifications

General Specifications:

Input: 90–240 VAC, 50–60 Hz, 2.3 Amp. Nom.

ETL and C-ETL Listed: Model ME 294 (9801427)

Domestic Model

Classification: Protective Class I Equipment

CE model Type BF Equipment

Enclosed equipment without protection against ingress of

water.

Equipment not suitable for use in the presence of a

flammable anesthetic mixture with air or with nitrogen oxide

Year 2000 Compliant Yes

Weight: 9.4 pounds

4,3 kg

Dimensions: 5 in (H) x 14.5 in (W) x 10 in (D)

12,7 cm (H) x 36,8 cm (W) x 25,4 cm (D)

Operating Temperature: $+50^{\circ}F$ to $+104^{\circ}F$

+10°C to +40°C

Humidity: Operating, 30% to 75% Relative Humidity at 104°F (40°C)

Nonoperating, 5 to 95% Relative Humidity, non-condensing

Storage Temperature: -40°F to 167°F

-40°C to 75°C

Storage Humidity: Storage, 30% to 90% Relative Humidity at 40° C, Non-

condensing

Storage Pressure: 700-900 mB

Environmental Disposal: The device contains lead in the form of solder used to

produce electrical contact between components. To avoid adverse environmental impact, utilize a disposal facility that performs complete incineration of the device at a temperature

in excess of 1000°C.

The shipping materials are fabricated of cardboard and may

be disposed of with other paper products.

Treatment timer:

Timer Accuracy: ± 0.5 minutes for times less than 5 minutes

 $\pm 10\%$ for times from 5 to 10 minutes

 ± 1.0 minute for times greater that 10 minutes

±5%, CE specification

Maximum Treatment Time: 60 minutes–electrical stimulation

Treatment Timer:

Treatment time counts down to zero when a time is set, or up to 60 or 30 minutes when no time is set. The digital timer indicates time in minutes and seconds. The timer also indicates the remaining or elapsed treatment time during the "Hold" period.

Waveform Specifications: Interferential Mode




Figure 3.8— Interferential Waveform

Waveform Type: Sinewave Polarity: None

Volts: 0–65 volts RMS, 1 Kohm load Current: 0–65 mA RMS, 1 Kohm load

Average current at maximum intensity

and frequency: 65 mA RMS

Maximum current density under 2"

diameter electrode. 3.2 mA/cm²

Frequency: Channel 1 = 4000 Hz

Channel 2 = 4000 to 4250

Hz variable frequency sine wave

Frequency Modulation: 1–15 Hz

80–150 Hz 1–150 Hz xx–xx Hz,

xx=any value from

1 to 250 Hz

Phase Duration: 125 μs

Available Amplitude

Modulation Options: Vector rotation

Available Channels: Channel pairs 1 & 2 or 3 & 4

Premodulated Mode



Figure 3.9—Premodulated Waveform

Waveform Type: Amplitude modulated sine wave

Polarity: None

Volts: 0–50 volts RMS, 1 Kohm load Current: 0–50 mA RMS, 1 Kohm load

Average current at maximum intensity

and frequency: 50 mA RMS

Maximum current density under

2" diameter electrode: 2.5 mA/cm²

Frequency: 4,000 Hz

Frequency Modulation: 1–15 Hz

80–150 Hz 1–150 Hz xx–xx Hz,

xx=*any value from*

1 to 250 Hz

Phase Duration: 125 µs internal sine wave

4-1,000 ms beat envelope

Available Amplitude

Modulation Options: Continuous

Surge

Reciprocation

Available Channels: All

Medium Frequency Mode




Figure 3.10—Medium Frequency (*Russian*) Waveform

Waveform Type: Burst modulated sine wave

Polarity: None

Volts: 0–50 volts RMS, 1 Kohm load Current: 0–50 mA RMS, 1 Kohm load

Average current at maximum intensity

and frequency: 50 mA RMS

Maximum current density under 2"

diameter electrode. 2.5 mA/cm²

Frequency: 2500 Hz, Burst at

10 ms on and 10 ms off

Frequency Modulation: No Phase Duration: 200 µs

Available Amplitude

Modulation Options: Continuous

Surge

Reciprocation

Available Channels: All

Biphasic (TNS) Mode




Figure 3.11—Biphasic (TNS)
Waveform

Waveform Type: Symmetrical biphasic square

wave

Polarity: None

Volts: 99 volts peak, 1 Kohm load

Current: 0 –99 mA peak, 1 Kohm load

Average current at maximum intensity

and frequency: 7.2 mA

Maximum current density under 2"

diameter electrode. 0.36 mA/cm² Frequency: 1–120 HzzHz

Frequency Modulation: No

Phase Duration: 50–300 µs

Available Amplitude

Modulation Options: Continuous

Surge

Reciprocation

Available Channels: All

High Volt Mode




Figure 3.12—High Volt Waveform

Waveform Type: Monophasic twin peak Polarity: Positive or negative

Volts: 500 volts peak, 1 Kohm load Current: 0–500 mA peak, 1 Kohm load

Average current at maximum intensity

and frequency: 1.2 mA at 120 pps with

1 Kohm load

Maximum current density under 2"

diameter electrode. 0.06 mA/cm² Frequency: 1–120 HzzHz

Frequency Modulation: No

Phase Duration: 8 µs at 50% Vmax Polarity: Positive or negative

Available Amplitude

Modulation Options: Continuous

Surge

Available Channels: Channel One only

Microcurrent Mode

•	Waveform Type:	Monophasic or biphasic square wave
or	Polarity:	Positive or negative or biphasic pulses
01	Volts:	1 Volt peak, 1 Kohm load
	Current:	10-990 μA peak, 1 Kohm load
Figure 3.13—Microcurrent Waveform	Average current at maximum intensity and frequency:	990 μΑ
	Maximum current density under 2"	
	diameter electrode.	$24.4 \mu\text{A/cm}^2$
	Frequency:	0.5-500 Hz
	Duty Cycle:	50%zHz
	Frequency Modulation:	No
	Pulse Duration:	1-1000 ms
	Available Amplitude Modulation Options:	Continuous
	Available Channels:	Channel Two only

Amplitude Modulation Specifications:

Vector rotation:	Interferential Mode Only	
	-50% amplitude modulation in anti phase with an eight second modulation period.	
Surge Mode:	Premodulated, Medium Frequency and Biphasic (TNS) Pulsed Modes	
Up ramp:	3 seconds	
Down ramp:	2 seconds	
Preset on/off times:	10 seconds on, 10 seconds off 10 seconds on, 20 seconds off 10 seconds on, 30 seconds off 10 seconds on, 40 seconds off 10 seconds on, 50 seconds off 10 seconds on, 60 seconds off	
Programmable On time:	1–240 seconds	
Programmable Off time:	1–240 seconds	
Reciprocation mode:	Premodulated, Medium Frequency and Biphasic (TNS) Pulsed Modes	
Up and down ramps: Reciprocation time:	1 second, <i>reciprocation only</i> 2–240 seconds, (On time = off time)	
Combine with Surge:	Use up and down ramps of surge program Use on/off times of surge program.	
Two timer option:	No	