



CARDIOTECH Standard vein probe

A medical device to visualize blood vessels in real-time

- Reveals the location and pattern of the blood vessels
- Non-invasive infrared detection technology makes it safe and hygienic

 $\sim h$

• Can be used with a broad range of patients; including infants, elderly and emergency patients

FEATURES

- Effective Enables to perform promptly vascular treat ments and procedures for the patients
- Accurate Enables to choose the right blood vessel based on real-time video of the blood vessel
- Precise Helps prevent repetitive and wrong injections
- Patient types Can be used with various patients such as infants and obese patients regardless of their skin tones

www.cardiacdirect.com 5963 Olivas Park Drive, Suite F Ventura, California 93003

Call our toll free number **1-888-354-2968**





Accessory




SLism-F: Specialized light source to see more clearly the veins of the face and other body parts

www.cardiacdirect.com 5963 Olivas Park Drive, Suite F Ventura, California 93003

Call our toll free number **1-888-354-2968**



Images



Details

| Dirobe Target | Subcutaneous blood vessels of various patients (including capillaries) |
|-------------------------------|--|
| - Probe Method | Non-invasive infrared technology |
| Probe Depth | Within 10 mm under the skin |
| Magnification Function | Up to 1 X, 2 X and 4 X |
| Automatic Adjustment Function | Adjusts the brightness automatically based on the targeted vein |
| Specialized Light Function | Shows the veins more clearly by using the specialized light source |
| Coperation Method | Easy to operate with touchscreen |
| Internal Battery Support | Continuous running time: 1,5 hours |

Application of Vein Probe



Intravenous injection

www.cardiacdirect.com 5963 Olivas Park Drive, Suite F Ventura, California 93003

Call our toll free number 1-888-354-2968