



ABPM-05 provides accurate information on blood pressure variability, overnight dipping and morning surge for reliable hypertension control.

Benefits

- meets BHS and AAMI accuracy criteria, recommended by dabl[®]
- proprietary CardioTech dual-valve system for excellent patient safety
- IP22 water-resistance rating
- sleep-friendly algorithm

manual programming option

7:29

ABPM-05

CardioTECH

- available in OEM
- free, all-inclusive software
- up to 48 hours of monitoring using 2 AA batteries

Cardio

normal adult cu

- up to 600 readings can be stored
- user friendly programming, with comprehensive analysis and reporting

Ambulatory blood pressure measurement gives information over and above conventional blood pressure measurement.

Why ABPM is needed?

- for more accurate BP profile
- to avoid the effects of white coat hypertension
- to record nightime BP, a major health risk indicator
- to lower the chance of misdiagnosis
- for better drug treatment

www.cardiacdirect.com 5963 Olivas Park Drive, Suite F Ventura, California <u>93003</u>

Call our toll free number **1-888-354-2968**



SPECIFICATIONS	ABPM-05	
dimensions	98*69*29 mm	
weight (excl. batteries)	190g	
memory	internal	
data transfer	USB optical cable	
power supply	2 AA batteries	
LCD display	40*22 mm	
numbers of keys	3 (start, event, day/night)	
operation	up to 51 hours	
measurement method	oscillometric, stepwise deflation	
measurement range	BP: 30-260 mmHg, pulse: 40-200 bpm	
passive accuracy	+/- 3 mmHg or 2% of the measured value	
sampling periods	4 (morning, active, passive, special	
storage capacity	600 measurements	
conformity certificates	CE, FDA	
validations	algorithm validated to BHS & AAMI protocols	

Free Software Solutions

Both for Busy GPs & Research Purposes

Both EasyABPM and CardioVisions programs display statistics for active, passive, morning and special periods. Statistics include systolic BP, diastolic BP, mean arterial pressure (MAP), pulse pressure (PP), pulse rate, minimum and maximum values, standard deviation (SD), percent time elevation (PTE), hypertension load (hyperbaric index), diurnal index (day/night index) and morning surge.

SPECIFICATIONS	EasyABPM	CardioVisions
device management	ABPM-05	ABPM-05
automated BP statistics	 Image: A second s	 Image: A second s
graphic & tabular display of data		✓
customizable sampling periods		 Image: A set of the set of the
printable report	 Image: A set of the set of the	 Image: A set of the set of the
pdf export		-
e-mail compatibility	-	 Image: A set of the set of the
password & login name	-	 Image: A set of the set of the
network use	-	✓