



Display

Type:	DPM7: 17" color TFT touchscreen DPM6: 12.1" color TFT touchscreen
Resolution:	DPM7: 1280 x 1024 pixels DPM6: 800 x 600 pixels
Waveforms:	DPM7: 12 selectable DPM6: 8 selectable

ECG (3 and 5-Lead)

Leads:	I, II, III, aVR, aVL, aVF, V
Gain Selection:	x0.125, x0.25, x0.5, x1, x2, auto
Sweep Speed:	12.5mm/sec, 25mm/sec, 50mm/sec
Bandwidth	Diagnostic Mode: 0.05-150Hz Monitor Mode: 0.5-40Hz Surgical Mode: 1-20Hz
Defibrillator Overload Protection:	Withstand 4000VAC/50Hz voltage in isolation against electrosurgical interference and defibrillation
Recovery Time:	<5sec
CMRR:	Diagnostic Mode: ≥ 90 dB Monitor Mode: ≥ 105 dB Surgical Mode: ≥ 105 dB (Notch filter set to off)

Heart Rate Meter

Measurement Range:	Adult: 15-300bpm Pediatric/Neonate: 15-350bpm
Accuracy:	± 1 bpm or $\pm 1\%$, whichever is greater
Resolution:	1bpm
Pacer Rejection:	When tested in accordance with the ANSI/AAMI EC13-2002: Sections 4.1.4.1 and 4.1.4.3, the heart rate meter rejects all pulses meeting the following conditions: Amplitude: ± 2 to ± 700 mV Width: 0.1 to 2ms Rise time: 10 to 100 μ s
Tall T-Wave Rejection:	When tested in accordance with the ANSI/AAMI EC13-2002 Section 4.1.2.1 c, the heart rate meter will reject all 100ms QRS complexes with less than 1.2mV of amplitude, and T-waves with T-wave interval of 180ms and those with Q-T interval of 350ms
Scaling Signal:	1mV $\pm 5\%$

ST Analysis

Adult/Pediatric Only	
Measurement Range:	-2.0mV to 2.0mV
Accuracy:	-0.8 to 0.8mV: ± 0.02 mV or $\pm 10\%$, whichever is greater
ST Adjust Scale:	60ms after J point, 80ms after J point (default: 60ms after J point)
ISO Adjust Scale:	4-200ms before R-Wave (default: 80ms) Step: 4ms
J Point Adjust Scale:	4-200ms after R-Wave (default: 48ms)

Pace Pulse

Pulse Indicator:	Pace pulses meeting the following conditions are marked by the PACE indicator: Amplitude: ± 2 to ± 700 mV Width: 0.1 to 2ms Rise time: 10 to 100 μ s
Pulse Rejection:	When tested in accordance with the ANSI/AAMI EC13-2002: Sections 4.1.4.1 and 4.1.4.3, the heart rate meter rejects all pulses meeting the following conditions: Amplitude: ± 2 to ± 700 mV Width: 0.1 to 2ms Rise time: 10 to 100 μ s

Arrhythmia Analysis

Adult/Pediatric Only

Asystole, ventricular fibrillation, ventricular tachycardia, pacer non-paced, pacer non-capture, ventricular rhythm, couplet, VT>2, bigeminy, trigeminy, R on T PVC, multiform PVC, irregular rhythm, missed beats, bradycardia, tachycardia

Respiration

Measurement Range:	Adult: 0-120bpm Pediatric/Neonate: 0-150bpm
Resolution:	1bpm
Accuracy:	7-150bpm: ± 2 bpm or $\pm 2\%$, whichever is greater 0-6bpm: undefined
Lead:	I or II (default: lead II)
Sweep Speed:	6.25mm/sec, 12.5mm/sec, 25mm/sec

Non-Invasive Blood Pressure

Measurement Method:	Oscillometric
Measurement Modes:	Manual, auto, stat
Connector Type:	Rectus
Units of Measure:	mmHg, kPa (user-selectable)
Resolution:	1mmHg
Systolic Range:	Adult: 40-270mmHg Pediatric: 40-200mmHg Neonate: 40-135mmHg
Diastolic Range:	Adult: 10-210mmHg Pediatric: 10-150mmHg Neonate: 10-100mmHg
Mean Range:	Adult: 20-230mmHg Pediatric: 20-165mmHg Neonate: 20-110mmHg
Accuracy:	Mean error is $\leq \pm 5$ mmHg Standard deviation is < 8 mmHg
Cuff Deflation Technique:	Step bleed



DPM™ | 6 | 7

Non-Invasive Blood Pressure (continued)

Initial Cuff Inflation:	Adult: 178 ±5mmHg Pediatric: 133 ±10mmHg Neonate: 87 ±5mmHg
Over Pressure Protection:	Double safety protection (hardware and software)
Pulse Rate Range:	40-240bpm
Pulse Rate Accuracy:	±3bpm or ±3%, whichever is greater

Invasive Blood Pressure

Measurement Range:	-50 to 300mmHg
Resolution:	1mmHg
Accuracy:	1mmHg or ±2%, whichever is greater
Zero Offset Range:	±200mmHg
Excitation:	5V DC, ±2%
Frequency Response:	DC to 12.5Hz ±1 Hz, -3db
Waveform Scales:	ART/Ao/UAP/BAP/FAP: 0 to 300mmHg PA: -6 to 120mmHg CVP/UVP: -10 to 40mmHg RAP/LAP/ICP: -10 to 40mmHg IBP1-IBP8: -50 to 300mmHg

Pulse Rate from Invasive Blood Pressure

Measurement Range:	25-350bpm
Resolution:	1bpm
Accuracy:	25-200bpm: ±1bpm or ±1%, whichever is greater 201-350bpm: ±2%

Pulse Oximetry

With Masimo SET® SpO₂

Measurement Range:	1-100%
Resolution:	1%
Accuracy:	±2% (70-100%, Adult/Pediatric, no motion) ±3% (70-100%, Neonate, no motion) ±3% (70-100%, Adult/Pediatric/Neonate, motion) 0-69% unspecified

Pulse Rate with Masimo SET® SpO₂

Measurement Range:	25-240bpm
Resolution:	1bpm
Accuracy:	±3bpm (no motion) ±5bpm (motion)

With Mindray SpO₂

Measurement Range:	0-100%
Resolution:	1%
Accuracy:	±2% (70-100%, Adult/Pediatric, no motion) ±3% (70-100%, Neonate, no motion) ±3% (70-100%, Adult/Pediatric/Neonate, motion) 0-69% unspecified

Pulse Oximetry (continued)

Pulse Rate with Mindray SpO₂

Measurement Range:	20-254bpm
Resolution:	1bpm
Accuracy:	±3bpm (no motion) ±5bpm (motion)

CO₂ with Mindray Sidestream

Measurement Range:	0-99mmHg
Resolution:	1mmHg
Accuracy:	0-40mmHg: ±2mmHg 41-76mmHg: ±5mmHg 77-99mmHg: ±10mmHg
Waveform Recognition:	0-40mmHg: ±2mmHg 41-76mmHg: ±5% reading 77-99mmHg: ±10% reading
Start-up Time:	<1min from start-up, module enters warming up status. 1min later, module enters ready-to-measure status (full accuracy mode)
Sampling Rate:	70ml/min or 100 ml/min (default: 100ml/min)
Auto-Zeroing Interval:	30sec, 10min and 30min after entering measurement mode and at every odd hour (1, 3, 5, 7, etc.) during operation after that
Respiration	
Measurement Range:	0-120bpm
Respiration Accuracy:	0-70bpm: ±2bpm >70bpm: ±5bpm

CO₂ with Oridion® Microstream®

Measurement Range:	0-99mmHg
Resolution:	Numeric: 1mmHg Waveform: 0.1mmHg
Accuracy:	0-38mmHg: ±2mmHg 39-99mmHg: ±5% + 0.08% × (reading - 38mmHg)
Waveform Recognition:	0-38mmHg: ±2 mmHg 39-99mmHg: ±5% of reading + 0.08% for every 1mmHg
Start-up Time:	30sec typical. Reaches 5% steady-state accuracy within 3min
Sampling Rate:	50ml/min: -7.5ml/min +15ml/min
Auto-Zeroing Interval:	At start-up, and every 12hrs thereafter
Respiration	
Measurement Range:	0-150bpm
Respiration Accuracy:	0-70bpm: ±1bpm 71-120bpm: ±2bpm 121-150bpm: ±3bpm

Anesthesia Gases

Sampling Rate:	Adult/pediatric: 120, 150, 200ml/min (user-selectable) (default: 120ml/min) Neonatal: 70, 90, 120ml/min (user-selectable) (default: 70/ml/min)
Sampling Delay Time:	<4sec
Refresh Rate:	1sec
Warm-up Time:	45sec to warm-up status 10min to ready-to-measure status
Normal Operating Conditions After Warm-up:	Ambient Temperature: 10 to 55°C (50 to 131°F) Ambient Pressure: 700–1200hPa Ambient Humidity: 10–95% RH, non-condensing
Measurement Range:	CO ₂ : 0-30% N ₂ O: 0-100% Des: 0-30% Sev: 0-30% Enf/Iso/Hal: 0-30% O ₂ : 0-100% AwRR: 2-100bpm
Resolution:	CO ₂ : 1mmHg AwRR: 1bpm
Accuracy:	CO ₂ : 0-1%: ±.1% 1-5%: ±.2% 5-7%: ±.3% 7-10%: ±.5% >10%: unspecified N ₂ O: 0-20%: ±2% 20-100%: ±3% Des: 0-1%: ±.15% 1-5%: ±.2% 5-10%: ±.4% 10-15%: ±.6% 15-18%: ±1% >18%: unspecified Sev: 0-1%: ±.15% 1-5%: ±.2% 5-8%: ±.4% >8%: unspecified Enf/Iso/Hal: 0-1%: ±.15% 1-5%: ±.2% >5%: unspecified O ₂ : 0-25%: ±1% 25-80%: ±2% 80-100%: ±3% AwRR: 2-60bpm: ±1bpm >60bpm: unspecified

Anesthesia Gases (continued)

Measurement Rise Time:	Sampling flow 120ml/min, using the DRYLINE™ water trap and neonatal DRYLINE™ 2.5m sampling line CO ₂ : ≤250ms N ₂ O: ≤250ms O ₂ : ≤600ms Hal/Iso/Sev/Des: ≤300ms Enf: ≤350ms
	Sampling flow 200ml/min, using the DRYLINE™ water trap and adult DRYLINE™ 2.5m sampling line CO ₂ : ≤250ms N ₂ O: ≤250ms O ₂ : ≤500ms Hal/Iso/Sev/Des: ≤300ms Enf: ≤350ms

Data Storage

Trend Data:	120hrs at 1min resolution 1hr at 1sec resolution
Alarm Events:	100 alarm events and associated waveforms (selectable waveform lengths: 8sec, 16sec or 32sec)
Arrhythmia Events:	100 arrhythmia events and associated waveforms (selectable waveform lengths: 8sec, 16sec or 32sec)
NIBP Measurements:	1,000 (systolic, diastolic, mean pressure, pulse rate and measurement time)
Full Disclosure Waveforms:	Maximum 24hrs (Specific storage time depends on the type and number of waveforms stored)

Recorder

Type:	Thermal array
Speed:	25mm/sec, 50mm/sec
# Traces:	3

Interfacing

Connectors:	1 AC power connector 1 RJ45 network connector, 100 BASE-TX (DPM6) 2 RJ45 network connector, 100 BASE-TX (DPM6) 4 USB 1.1 connectors (DPM6) 10 USB 1.1 connectors (DPM7) 1 nonstandard USB SMR connector 1 50-pin CF revision 2.0 connector 1 standard DVI-D video interface connector 1 BNC connector 1 equipotential grounding connector 1 RJ11 defib sync connector 1 CIS connector (DPM6) 1 Micro D connector (DPM6)
-------------	---



Battery

Type:	Rechargeable lithium ion
Number of Batteries:	2
Run Time:	Lithium ion: 2hrs using 2 new, fully charged batteries and monitoring ECG, SpO ₂ and auto NIBP measurements every 15min at 25°C
Recharge Time:	6hrs maximum

Physical Dimensions

DPM7 Monitor Size:	37cm(H) x 40cm(W) x 19.3cm(D) 14.6"(H) x 15.7"(W) x 7.6"(D)
DPM6 Monitor Size:	29.7cm(H) x 33.6cm(W) x 18.6cm(D) 11.7"(H) x 13.2"(W) x 7.3"(D)
DPM7 Monitor Weight:	Less than 14.5kg (31.9lbs) including 17" touchscreen display, MPM, AG module, 2 lithium batteries, recorder, and CF components
DPM6 Monitor Weight:	7.2kg (15.8lbs) including 12.1" touchscreen display, MPM, and ECG, NIBP, SpO ₂ accessories

Environmental

Operating Temperature:	0°C to 40°C (DPM6 main unit/MPM/IBP module/recorder) 5°C to 40°C (DPM7 main unit/MPM/IBP module/recorder) 0°C to 40°C (Microstream CO ₂ module) 5°C to 35°C (Sidestream CO ₂ module) 10°C to 40°C (AG module)
Storage Temperature:	-20°C to 60°C (main unit/MPM/IBP module/recorder) -20°C to 60°C (Microstream CO ₂ module) -20°C to 60°C (Sidestream CO ₂ module) -20°C to 70°C (AG module)
Operating Humidity:	15% to 95%, non-condensing (main unit/MPM/IBP module/recorder/Microstream CO ₂ module/Sidestream CO ₂ module/AG module)
Storage Humidity:	10% to 95%, non-condensing (main unit/MPM/IBP module/recorder/Microstream CO ₂ module/AG module) 15% to 95%, non-condensing (Sidestream CO ₂ module)
Operating Atmospheric Pressure:	425-809mmHg (main unit/MPM/IBP module/recorder) 430-795mmHg (Microstream CO ₂ module) 428-790mmHg (Sidestream CO ₂ module) 525-900mmHg (AG module)
Storage Atmospheric Pressure:	120-809mmHg (main unit/MPM/IBP module/recorder) 430-795mmHg (Microstream CO ₂ module) 428-790mmHg (Sidestream CO ₂ module) 525-900mmHg (AG module)

Power Requirements

AC Voltage:	100 to 240VAC, 50/60Hz
DPM7 Current:	2.8 to 1.6A
DPM6 Current:	2.5 to 1.4A

Safety

Type of Protection:	Class I (main unit and secondary display only)
Degree of Protection:	Sidestream CO ₂ module/Microstream CO ₂ module/AG module: BF MPM/IBP module: CF
Protection Against Ingress of Fluids:	Not protected

Mindray DS USA, Inc.

800 MacArthur Blvd., Mahwah, NJ 07430

Tel: 1.800.288.2121 Tel: 201.995.8000 Fax: 1.800.926.4275

www.mindray.com

mindray™
NORTH AMERICA