

Medical

a subsidiary of
VIASYS Healthcare

The precision Peak Flow
Meter with an immaculate
pedigree



With EN13826
European Scale

**Features**

- *Designed specifically to meet EN13826 scale standard*
- *Movable asthma management colour zones that won't wash off*
- *Large, high visibility scale*
- *In-built internal check valve, avoiding cross infection concerns*
- *Meets and exceeds all international performance standards*
- *Customisation possible for promotional purposes*
- *Available with either European or ATS scale*

Focus on the Future




Medical

25 years experience of design, testing and improving the world's leading Peak Flow Meters (PFM) has led to the new generation MicroPeak. Designed to meet the new European (EN13826) standards in performance and quality including a high visibility scale and in-built Colour-Zone asthma management system that won't wash off (see illustration below).

Unlike other PFM's the advanced design of MicroPeak means dismantling is avoided, thus protecting the integrity of the internal mechanism, also an internal check valve avoids cross infection concerns.

MicroPeak comes complete with instructions and a 4 week recording chart, and can be offered fully customised for promotional purposes too.

When accuracy is important in clinical asthma management, then compromise is not an option, choose MicroPeak the PFM with an immaculate pedigree.



(Colour-Zone asthma management system)

MicroPeak EU scale ▶

Technical Specifications

Range: 60-900 l/min **(EU, (EN13826) (Cat. No. MPE8200EU))**
60-900 l/min **(ATS, (Cat. No. MPE7200))**

Scale Increment: 10 l/min

Accuracy: <+/- 5% or 10 l/min

Repeatability: <+/- 2%

Calibration: 100% (individual unit)

Frequency Response: <100Hz @ 350 l/min


Dynamic State: Critically damped @ 350 l/min

Standard Compliance: ERS Statement 1997
(Eur Respir J 1997;9 Suppl.24)
AS/NZS 4237:1994
NAEP Guidelines 1991

Integral Check Value Expiration only



▲ MicroPeak ATS scale



Peak flow record chart showing recovery from acute asthma with pronounced morning dips.

These Peak Flow Meters are part of an extensive range of respiratory measurement devices manufactured by Micro Medical Ltd and are offered as Catalogue No MPE8200EU MicroPeak European scale or Catalogue No MPE7200 MicroPeak ATS (American Thoracic Society) scale.

Micro Medical Ltd pursues a policy of continuing improvement in design, production and performance of it's products, the right is therefore reserved to vary details at any time without notice.



Micro Medical Limited
PO BOX 6, Rochester, Kent, ME1 2AZ, UK
Telephone 01634 893500
Fax 01634 893600
International +44 1634 893500
Email micromedical@viasyshc.com
www.micromedical.co.uk



Local Representative:

