



Propaq LT | Highly Deployable Continuous Monitor

Introducing the Innovative New  
**Welch Allyn Propaq<sup>®</sup> LT**

It won't just change how you monitor patients.  
It'll change where you do it.



**WelchAllyn<sup>®</sup>**

Advancing Frontline Care<sup>™</sup>

# The New Welch Allyn Propaq® LT

Monitor more patients in more places than ever before.



Reliability, durability and affordability combine to give you perhaps the most versatile monitor on today's market.

The new Propaq LT is the perfect monitor for today's fast-paced frontline hospital environment—letting you extend continuous monitoring to areas in the hospital where technical complexity or cost made it impossible in the past. Because it's so easy to use and configure, the new Propaq LT is the perfect choice for nurses on general care floors. And with legendary Propaq quality, you know it's rugged and reliable enough to use in the ED and for interdepartmental transport.

With wireless connectivity and spot-check mode options, the new Propaq LT is like getting three monitors—ambulatory, bedside and transport—for the price of one.




The one monitor you can deploy when—and where—you need it most.



Equipped with the optional Propaq LT large color display and charging cradle, the Propaq LT is ideal for use as a bedside or procedure monitor. And with a revolutionary new design that's light enough to fit on an IV pole, rugged enough for air transport, and versatile enough for pediatric applications, the Propaq LT gives you the power to monitor more patients in more places than ever before.



SHOWN ACTUAL SIZE:  
The Welch Allyn Propaq LT  
weighs about two pounds.



Tab/Wave	80	127/71	12	97
Time	HR/min	NIBP mmHg (86)	Resp/min	SpO2
12:30	80	n/a	12	97
12:29	80	n/a	12	97
12:28	80	n/a	12	97
12:27	80	n/a	12	97
12:26	80	127/71 (86)	12	97
12:25	80	n/a	12	97
12:24	80	122/73 (86)	12	97

## Key Features:

- Shock resistance of up to 75 g—enough to withstand a 6-foot drop
- Heart rate, 3- or 5-lead ECG, SpO<sub>2</sub>, NIBP, respiration
- Neonatal, pediatric and adult patient modes
- Continuous monitoring and spot-check modes
- 24-hour tabular and graphical trending
- Storage of up to 20 full-disclosure data snapshots
- Wireless connectivity option to Welch Allyn Acuity® Central Station via FlexNet™
- Designed for use in large- and small-body fixed- and rotary-wing aircraft

Quickly toggle between three display modes that are right for your care area (five to choose from). Shown are two-waveform, tabular trend with waveform, and large numeric (pictured above) displays.



# Designed to meet the needs of today's hospital environments.



Large color display shown on wall mount.

## Large Color Display Option

The Propaq LT conveniently connects to a large color display with the choice of 4- or 9-waveform display modes, making it perfect for bedside or procedure applications. And with tabletop and wall mounting options, you have the flexibility to use your large color display where you need it most. Off-the-shelf displays may be purchased separately for use with this option.



Choose from 4- or 9-waveform display modes and easily see vital signs numerics, waveforms and alarm limits.



## Incorporating the Latest Technologies

The Propaq LT can communicate with any Windows 2000/XP PC via USB port, allowing the user to custom-configure the monitor and print tabular trends and snapshots of patient data. And with a wireless connectivity option to the Acuity Central Station, the Propaq LT can integrate into your hospital's existing FlexNet network.



All Propaq vital signs monitors come standard with Smartcuf technology\*, a new standard in NIBP accuracy. This patented software that improves NIBP performance and accuracy in the presence of motion artifact.

*\*SmartCuf available in future Propaq LT software releases coming in 2006.*



Only Nellcor® pulse oximetry sensors, including *OxiMax*® patient-dedicated adhesive sensors, should be used with the Nellcor pulse oximetry option. Nellcor is a registered trademark of Nellcor Puritan Bennett Inc.

Welch Allyn  
8500 S.W. Creekside Place  
Beaverton, OR 97008-7107  
(t) 503.530.7500 (f) 503.526.4200  
Toll-Free: 800.289.2500  
Technical Support: 800.289.2501  
www.welchallyn.com



Small enough to fit on a nightstand, light enough to carry around, simply drop the compact new Propaq LT into the self-charging bed rail cradle for use in bedside applications.

Contact your local Welch Allyn sales representative for more information today.

**WelchAllyn®**

Advancing Frontline Care™

© 2006 Welch Allyn

All rights reserved. Acuity, Propaq and Smartcuf are registered trademarks and FlexNet is a trademark of Welch Allyn. Other product and company names mentioned herein are for identification purposes only and may be the trademarks or registered trademarks of their respective companies.

Printed in USA SM2832 Rev D