

COMEN

Share With The World.



C50

Multi-parameter Patient Monitor



CE 0434

C50 Multi-parameter Patient Monitor Features

- 10.4" TFT Touch Screen.
- Operating methods: Trim knob/ keybutton, touch screen.
- Advanced patient file management system,
- High capacity storage, power-off data protection.
- Inclined screen for easy observation.
- **Standard Configuration:**
5 lead ECG, HR, COMEN SpO₂, PR, NIBP, RESP, Dual TEMP
- **Optional Configuration:**
Nellcor SpO₂, Maximo SpO₂, Dual IBP, EtCO₂, Thermal printer



- Fixed handle, easy to move



- 4000mAh Li-ion batteries, above 4 hours working time



- Auto dust shield to prevent from foreign matters



- Handwriting pen with slot Interfaces: SD card, USB, Network, VGA.



- Mounting solutions: vertical hanging, wall mounting, rolling stand.



- Complying with CE standards, integrated modular design and aluminum alloy structure gives better guarantee to product safety & stability.
- International pop design, splendid appearance, best craftsmanship.



- IPX1 waterproof design, good performance in various environment.



- I-Klok™ Intelligent alarm management, automatic-identification of alarm level. Self-adjust proper alarm time to reduce false alarms.

- ARR analysis, RESP Apnea Alarm, OxyCRG.

- CMRR ≥ 105 dB, outstanding ECG anti-interference capability.



C80 ICU Monitor



- CardioTec™ 12-lead ECG/HR, OxiMax® SpO₂, AcuTec™ NIBP, RESP, TEMP, 12.1" TFT touch screen, handwriting pen, EtCO₂ socket, Li-ion battery.
- Option: 2-IBP, EtCO₂, AG, ICG (CO) , printer, pediatric Neonatel accessories.
- IPX1 level waterproof design.
- I-Klok™ intelligent alarm management for all the parameters, auto-identification for alarm level.

C60 NICU Monitor



- According to the feature of neonatal HR, we use the latest ExNeo™ technology to guarantee the accuracy of measurement.
- World leading SpO₂ technology, OxiMax® SpO₂.
- SatSeconds™ intelligent alarm management technology to reduce false alarms, accurate patient information transfer.
- According to the low blood and low perfusion feature of neonate, we use Adap-DSP™ system to provide more accurate measurement.
- Unique neonate measurement mode, soft & hard over-pressure protection.
- Bio-compatibility approved accessories for neonate.

Shenzhen Comen Medical Instruments Co., Ltd.

Add: Floor 7, Block 5, 4th Industrial Park of Nanyou, Nanshan District, Shenzhen City 518052, China

Tel: +86-755-2640 8879 2641 9446 Fax: +86-755-2643 1232

Website: www.szcomen.com E-mail: info@szcomen.com