

Vetus EQEquine Ultrasound Image System

Efficiency & Quality









Q | Search Mindray Animal Medical

www.mindrayanimal.com service@mindrayanimal.com







Mindray Animal Medical has been developing innovative and diversified ultrasound veterinary solution since its inception. For years, Mindray Animal Medical always strives to be supportive of healthcare of varied species, including companion animals, rare animals, farm animals and any other ones.

Vetus EQ is a portable and practical color doppler ultrasound system customized for the equine diagnosis. Based on the revolutionary technology of ZST⁺ platform, Vetus EQ provides the professional solutions for ultrasonic examination of horses.

Concise & Practical

Innovative Design

- Intuitive touch interface
- Sealed touchpad
- Seamless design control panel

Sturdy Construction

- Touchpad design prevent animal hairs
- Resist corrosion from harsh cleaning agents




 15.6-inch high revolution color LED monitor, narrow, border, wide view

 12.3-inch high sensitivity anti-glare color touch screen



Advanced ZST⁺ Imaging Platform

ZST⁺ has three core technologies: Zone Imaging, Zone Focusing, and Zone Processing, which greatly improved the imaging level of an ultrasound system in all directions, help equine practitioners see more precisely, comprehensively and quickly.




Dedicated Equine Applications



iWorks dedicated procedure

Standard and simplified the workflow dedicated for equine application reduce 50% exam time, 80% keystrokes



MSK Application

L13-3Ns

Wide frequency range linear, for equine tendon application



Abdomen Application

C5-1s

Convex array transducer for equine abdominal scanning in big animals






Reproduction Application

6LE5Vs

Intrarectal linear array transducer




Cardiology Application



P4-2s
Excellent penetration
phased array for equine
cardiology application





Portable & Durable

iCover—Protective Cover Set

Prevent dust, dirt and splash, anti accidental fall, pump and smash.



U-bank

Main unit supports up to 8 hours of continuous scanning with U-bank.



Easy Backpack

Suitable and sturdy backpack for easy carriage especially in outside work conditions. Dedicated design of inner pockets for transducers, adapter modules, coupling gels and other tools.



Dedicated & Professional

Tendon Application



iScape

iScape provides wide-view image to help measure the equine tendon length and evaluate its status.



iNeedle+

iNeedle⁺ provide the equine practitioner with an easy tool for the operations: injection, aspiration, or biopsy.

Specific and Professional Presets




Focused on the equine ultrasound diagnosis demand, Vetus EQ developed a series of professional preset for different application scenarios:

- Equine MSK
- Equine Abdomen
- Equine Cardiac
- Equine Reproduction


All above increase the efficiency and quality of equine pracitioners' work.



Equine Tendon



Equine Cardiac TVD




Equine Proximal Metacarpal Region



Tendon Calcification



Equine Tendon Hydrops



Equine Lumbar Muscle



Equine Abdominal Aorta Flow



Equine Spine